


Welch Allyn® Spot® Vision Screener Model VS100



An easy to use autorefractor that accurately identifies refractive errors and ocular misalignment in children and adults.



Designed for today's vision screening environment

Three distinct screening modes fit your club preferred workflow and data management needs



Instant Screening



Select from Queue



Enter Subject Data & Screen




Use Instant Screening Mode

Setup

- Ensure that room lighting is appropriate for vision screening. Turn off lights and close blinds or curtains, if possible. The child's back should be positioned to the interfering light source.
- Avoid direct light into the patient's eyes from: ambient light sources, exterior lights, incandescent and halogen lights.
- Explain the process to child and parent/guardian
- Ergonomic design is comfortable for screener while capturing from a non-invasive 3.3 feet distance

Proper Setup



Alternate Operator Position

Operator seated in a chair with device placed on top of tripod platform. Slowly move device forward and back on platform. Operator assisted by helper.



Capturing Results

- From the Home Screen, select the appropriate screening mode. *Note: Screen will dim after 60 seconds, just tap screen to awake.*
- Slowly rotate the device upward to meet both of the subject's eyes and lean forward and backward keeping the pinwheel steady.
- Position yourself with one foot in front of the other, slowly rock forward and back until the screen turns grey, indicating you are in the capture range.
- Alternately, position yourself in a chair with device placed on top of tripod platform. Slowly move device forward and back on platform until the screen turns grey, indicating you are in the capture range.
- Auto capture feature insures screening conditions are optimal. It ONLY captures when the patient is looking at the device and the user has the right conditions (distance, lighting, etc.).



On-screen instructions guide the user with a "too close" or "too far" indicator


Capturing Results

Keep the vision screener device steady until the screening wheel appears, indicating the capture process is underway. Capture is complete when results screen is shown.



Managing Results

Immediate on-screen results offers either a “Complete Eye Exam Recommended” or “All Measurements in Range” so no interpretation of data is required. Also lists issues detected if referred.



When results equal or exceed the criteria threshold settings, the results are:

Out-of-Range/Refer =
“Complete Eye Exam Recommended”

If below the criteria results are:

In-Range/Pass =
“Screening Complete all measurements in range”